

# H - 80 Series

	H-80 AS	H-80 A	H-80 M
Adjustment method of oxygen concentration	Auto	Auto	Manual
Flow	Monitorable, adjustable	Monitorable, adjustable	Monitorable, adjustable
Flow range	2L-80L/min	2L-80L/min	2L-80L/min
Oxygen concentration adjustment range	100% 21%-	100% 21%-	100% 21%-
• Temperature output range	29° C – 37° C	29° C – 37° C	29° C – 37° C
• Temperature adjustment gear	9 steps adjustable	9 steps adjustable	9 steps adjustable
Disinfection free conception	Yes	Yes	Yes
Screen size	3-5 inches	3-5 inches	3-5 inches
Trend review function	1day, 3days, 7days	1day, 3days, 7days	1day, 3days, 7days
Oxygen concentration	Monitorable, adjustable	Monitorable, adjustable	Monitorable
Temperature	Monitorable, adjustable	Monitorable, adjustable	Monitorable, adjustable
• sPEEP	Monitorable	Monitorable	Monitorable
• Respiratory rate	Monitorable	Monitorable	Monitorable
Automatic reminder of filter replacement	Yes, settable	Yes, settable	Yes, settable
Automatic water filing function	Yes	Yes	Yes
• Auto flow	Yes	/	/
• Smart flow	Yes	/	/
Humidity compensation	3 +-, +3 step adjustable 7	3 +-, +3 step adjustable 7	3 +-, +3 step adjustable 7
• sPEEP	steps adjustable 0-4	/	/
• Warm stand by	Yes	Yes	Yes



## H-80 Series

High flow nasal cannula (HFNC) oxygen therapy device



2A El-Kaliefa El-Mamon St.  
Masr Elgededa , Cairo ARE  
Tel & Fax : 0224537297  
[www.iand-a.com](http://www.iand-a.com)



KM 25 Alexandria- Cairo Desert Rd  
Amriya , Alexandria , Egypt  
Tel : +20 (3) 4700083 / 9613699  
Fax : +20 (3) 4703312  
[www.Pharco.org](http://www.Pharco.org)  
[Pharco@pharco.com.eg](mailto:Pharco@pharco.com.eg)

